



BioSurface Technologies Corporation

Microscopy and Image Analysis Flow Cells

[All BioSurface Technologies Flow Cells are Completely Autoclavable and Reuseable]

Flow Cell Accessories

DESCRIPTION

PRICING



FC 33 Bubble Trap Single Channel

\$75.00*

Collects gas bubbles from the influent supply line prior to the Flow Cell
Barbed fittings accommodate size 16 (1/8"=3.2 mm ID) tubing.
Autoclavable and re-useable.



FC 34 Bubble Trap 4-channel

\$195.00*

Collects gas bubbles from the influent supply line prior to the Flow Cell
Barbed fittings accommodate size 16 (1/8"=3.2 mm ID) tubing.
Autoclavable and re-useable.



FB 50 Glass Flow Break (single channel)

\$36.00*

Barbed fittings accommodate size 16 (1/8"=3.2 mm ID) tubing.
Autoclavable and re-useable.